


FINNPIPETTE DIGITAL

INSTRUCTIONS FOR USE



Thermo

ELECTRON CORPORATION

Thermo Electron Oy
P.O.Box 100, Fin-01621 Vantaa, Finland
Tel. +358-9-329 100, fax -358-9-3291 0414
www.thermo.com/finnpipette

PIPETTE DESCRIPTION

The FINNPIPETTE DIGITAL is a continuously adjustable, general purpose micropipette for sampling and dispensing accurate amounts of liquids. It operates on air displacement principle (= air interface) and uses detachable, disposable tips.

The adjusted delivery volume is displayed in clear digits on a readout window in the handle.

FINNPIPETTE DIGITAL pipettes cover a volume range of 0,5-10 ml.

Code	Volume range	Scale division
4027000	0,5 - 10 μ l	0,1 μ l
4027250	2 -20 μ l	0,1 μ l
4027010	5 - 40 μ l	0,5 μ l
4027260	20-200 μ l	1,0 μ l
4027020	40 - 200 μ l	1,0 μ l
4027030	200 - 1000 μ l	5,0 μ l
4027040	1 - 5 ml	50 μ l
4027050	2-10 ml	100 μ l

All models are equipped with a built-in tip ejector.

RAW MATERIALS

The FINNPIPETTE DIGITAL is made of mechanically durable and chemically resistant materials.

DIGITAL DISPLAY (Fig.1)

The delivery volume of the pipette is displayed on read-out window in the handle.

DESCRIPTION OF TIPS (Page16)

Finntips are recommended for use with the FINNPIPETTE DIGITAL.

They are made of natural coloured polypropylene, generally regarded as the only contamination free material suitable for tips.

Finntips are also autoclavable (121°C).

PIPETTE OPERATION

SETTING THE DELIVERY VOLUME (Fig.2)

To increase the delivery volume turn the operating button counterclockwise and to decrease turn clockwise (operating button facing you) .

Make sure the desired delivery volume clicks into place and the digits are completely visible.

Caution: Accuracy and precision are valid only within the nominal volume range.

Do not set volumes outside the range. Using excessive force to turn the operating button outside the range may jam the mechanism and eventually cause defects.

TIP EJECTION (Fig.3)

Direct the pipette towards a suitable waste receptacle.

Press the tip ejector down and the tip will be cleanly ejected.

PIPETTE TECHNIQUES

FIGURES 4-7:

- A = READY POSITION**
- B = FIRST STOP**
- C = SECOND STOP**

The operation is controlled by the operating button.

For best possible accuracy:

Operate the operating button slowly and with steady speed at all times, especially when working with liquids of high viscosity.

Never let the operating button snap back. Make sure the Finntip is firmly attached to the tip cone of the pipette and that there are no foreign bodies inside the tip itself.

Wet a newly attached tip with the solution being pipetted before any actual pipetting takes place.

This is done by filling and emptying the tip 2-3 times.

Hold the pipette vertically during intake of the solution.

The temperature of the tip and pipette should be equalized to that of the solution.

For maximum hand comfort hold the pipette lightly in your hand with the grippy resting over your index finger.

FORWARD TECHNIQUE (Fig.4)

1. Depress the operating button to the first stop.
2. Immerse the tip slightly about 1 cm below the surface of the solution and release the operating button slowly. The tip has now been filled. Withdraw the tip from the solution. Wipe any drops from the outside of the tip without touching the orifice.
3. Deliver the liquid by gently depressing the operating button to the first stop. Then after a delay of about a second continue to depress the operating button all the way down to the second stop. The tip has now been emptied. Drop formation is avoided by touching the side of the vessel with the tip. Withdraw the tip.
4. Release the operating button to the ready position.

If necessary change the tip and continue with your pipetting.

REVERSE TECHNIQUE (Fig.5)

The reverse technique is suitable for pipetting solutions with high viscosity and/or solutions that tend to foam easily. This technique is also recommended for pipetting very small volumes.

1. Depress the operating button all the way to the second stop.
2. Immerse the tip slightly about 1 cm below the surface of the solution and release the push button slowly. The tip has now been filled. Withdraw the tip from the solution with the tip sliding along the vessel wall.
3. Deliver the liquid by gently depressing the operating button to the first stop. The tip has now been emptied. Drop formation is avoided by touching the side of the vessel with the tip. Withdraw the tip. Some liquid will remain inside the tip and should not be included in the delivery.
4. The remaining liquid is either discarded with the tip or pipetted back into its original container.

Please note: if bigger volumes than 4 ml are pipetted with reverse technique there is no room in the tip for the extra remaining liquid. Normally reverse technique is not used for big volumes like this.

However bigger volumes than 4 ml of foaming liquid can easily be pipetted as follows: Use forward technique steps 1 and 2 to fill the tip.

Deliver by gently depressing the operating button to the first stop. Drops outside the tip should be included in the delivery.

REPETITIVE TECHNIQUE (Fig.6)

The repetitive technique offers a rapid and simple procedure for repeated deliveries of the same liquid with the same volume.

1. Depress the operating button all the way down to the second stop.

2. Immerse the tip slightly about 1 cm below the surface of the solution and release the push button slowly.

The tip has now been filled.

Withdraw the tip from the solution with the tip sliding along the vessel wall to avoid drops.

3. Deliver the liquid by gently depressing the operating button to the first stop.

The desired volume has not been delivered.

Hold the operating button at the first stop.

Some liquid will remain inside the tip and should not be included in the delivery. Drops outside the tip should be included in the delivery.

4. Immerse the tip slightly below the surface of the solution and release the operating button slowly

The tip has been refilled.

5. Continue repeating procedures 3 and 4.

PIPETTING WHOLE BLOOD (Fig.7)

e.g. deproteinization in blood glucose determination)

Use forward technique procedures 1 and 2 to fill the tip with blood.

Wipe the tip carefully with a dry and clean tissue.

1. Immerse the tip into the reagent and depress the operating button to the first stop, making sure the tip is well below the surface.

2. Release the operating button slowly to the ready position.

The tip has now been filled with reagent.

Do not lift the tip out of the solution.

3. Depress the operating button to the first stop and release slowly.

Keep repeating this procedure until the interior wall of the tip is clear.

4. Finally, depress the operating button all the way to the second stop to completely empty the tip.

MAINTENANCE

When the Finnpiquette Digital is not in use make sure it is safely stored in a vertical position.

SHORT TERM CHECKING

The pipette should be checked at the beginning of each day for dust and dirt on the outside surfaces of the pipette. Particular attention should be paid to the tip cone of the pipette.

Solvents should not be used for cleaning the pipette.

LONG TERM MAINTENANCE (0,5-5 ml) (Fig.8)

The piston and the cylinder should be checked at least twice a year if the pipette is in daily use.

Opening the pipette

1. Depress the tip ejector pusher.
2. Insert the tooth of the opener into the opening at the base of the tip ejector pusher.
3. Pull away the tip ejector shaft and tip ejector pusher.
4. Remove the tip cone by turning it in a counterclockwise direction with the service tool provided in the package.
5. Pull out the piston.
6. Remove the O-ring from the tip cone. In model 5-40 μ l the two O-rings are located deep inside the tip cone. They can be removed by using the thin end of the piston.
7. Clean the piston, the piston spring and the O-ring with a dry, napless cloth.
8. Check the cylinder for foreign particles. Please note that the cylinder must not be greased.
9. Grease the cleaned parts with high vacuum stop cock lubricant.
10. Reassemble the parts.

LONG TERM MAINTENANCE (10ml) (Fig.9-10)

If the pipette is used daily it should be checked at least twice a year. The servicing procedure starts with disassembly of the pipet.

1. First remove the Tip Ejector pusher by pushing it all the way down and lifting the base with a screwdriver blade or with the tooth of the Service Tool to release the snap joint.
2. Next remove the Tip Ejector part 2 from part 1 using the maintenance pliers to release the snap joint.

3. Remove the Cylinder by pressing the tip ejector part 1 firmly towards the cylinder. This action releases the snap joint and you can pull the cylinder away.
4. Clean the O-ring and cylinder. Regrease the O-ring and the cylinder.
5. Assemble the parts in the opposite order to removal. All joints are a snap fit and are made simply by pushing the parts together by hand. Be careful not to bend the pipette during assembly because this could damage the snap joints or the piston.

CAUTION!

The Finnpiquette is designed to allow easy in-lab service. If you, however, want to send the pipette to us or to our local representative for service, please enclose a list of any infectious, radioactive or otherwise hazardous materials that have been pipetted. Also please note that the postal authorities in your country may limit the sending of contaminated material by mail.

CALIBRATION

All Finnpipettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. Normally, the pipettes do not need adjustment, but they are constructed to permit recalibration and adjustment for liquids of different temperature and viscosity.

DEVICE REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range	readable graduation
under 10 μl	0.001 mg
10-100 μl	0.01 mg
above 100 μl	0.1 mg

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ($\pm 0.5^\circ\text{C}$) temperature of water, pipette and air between 20°C to 25°C .

The relative humidity must be above 55%. Especially with volumes under 50 μl the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

CHECKING THE CALIBRATION

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume or 10% of maximum volume, whichever is higher. E.g. Finnpiquette 0.5-10 μl is tested at 10 μl and 1 μl . A new tip is first pre-wetted 3-5 times and a series of ten pipettings are done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume. Measuring volumes taken from balance is not allowed. If the calculated results are in the limits, the calibration of the pipette is correct.

Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the accuracy (A) and precision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the results are in the limits of Table 1, then the calibration of the pipette is correct. Otherwise the pipette must be adjusted and checked again.

Range	Volume μl	Accuracy		Precision	
		μl	%	s.d.* μl	cv%
0,5-10 μl	10	± 0.10	± 1.0	0.08	0.8
	1	± 0.04	± 3.5	0.03	3.0
2-20 μl	20	± 0.200	± 1.0	0.080	0.4
	2	± 0.060	± 3.0	0.030	1.5
5-40 μl	40	± 0.240	± 0.6	0.120	0.3
	5	± 0.125	± 2.5	0.100	2.0
20-200 μl	200	± 1.2	± 0.6	0.6	0.3
	20	± 0.6	± 3.0	0.3	1.5
100-1000 μl	1000	± 5.0	± 0.5	2.0	0.2
	100	± 1.5	± 1.5	0.6	0.6
1-5 ml	5000	± 25.0	± 0.5	10.0	0.2
	1000	± 15.0	± 1.5	5.0	0.5
2-10 ml	10000	± 50.0	± 0.5	20.0	0.2
	2000	± 20.0	± 1.0	6.0	0.3

ADJUSTMENT (Fig.10)

Adjustment is done with the service tool.

1. Place the service tool into the openings of the calibration nut at the top of the handle.
2. Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.
3. After adjustment check the calibration according to the instructions above.

FORMULAS FOR CALCULATING RESULTS

Conversion of mass to volume

$$V = (w + e) \times Z$$

V = volume (μl)

w = weight (mg)

e = evaporation loss (mg)

Z = conversion factor for mg/ μl conversion

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6 mg = 0.2 mg/s). Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2 mg (10s x 0.2mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary. The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032 $\mu\text{l}/\text{mg}$ at 22°C and 95 kPa. See the conversion table on page 10.

Accuracy (systematic error)

Accuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0$$

A = accuracy
 \bar{V} = mean volume
 V_0 = nominal volume

Accuracy can be expressed as a relative value:

$$A\% = 100\% \times A / V_0$$

Precision (random error)

Precision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$s = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

s = standards deviation

v = mean volume

n = number of measurements

cv is the relative value of standard deviation.

$$cv = 100\% \times s / \bar{v}$$

CONVERSION TABLE

Value of the conversion factor Z ($\mu\text{l}/\text{mg}$), as a function of temperature and pressure, for distilled water.

Temperature °C	Air pressure hPA (mbar)					
	800	853	907	960	1013	1067
15	1.0018	1.0018	1.0019	1.0019	1.0020	1.0020
15.5	1.0018	1.0018	1.0019	1.0020	1.0020	1.0021
16	1.0019	1.0020	1.0020	1.0021	1.0021	1.0022
16.5	1.0020	1.0020	1.0021	1.0022	1.0022	1.0023
17	1.0021	1.0021	1.0022	1.0022	1.0023	1.0023
17.5	1.0022	1.0022	1.0023	1.0023	1.0024	1.0024
18	1.0022	1.0023	1.0024	1.0024	1.0025	1.0025
18.5	1.0023	1.0024	1.0025	1.0025	1.0026	1.0026
19	1.0024	1.0025	1.0025	1.0026	1.0027	1.0027
19.5	1.0025	1.0026	1.0026	1.0027	1.0028	1.0028
20	1.0026	1.0027	1.0027	1.0028	1.0029	1.0029
20.5	1.0027	1.0028	1.0028	1.0029	1.0030	1.0030
21	1.0028	1.0029	1.0030	1.0030	1.0031	1.0031
21.5	1.0030	1.0030	1.0031	1.0031	1.0032	1.0032
22	1.0031	1.0031	1.0032	1.0032	1.0033	1.0033
22.5	1.0032	1.0032	1.0033	1.0033	1.0034	1.0035
23	1.0033	1.0033	1.0034	1.0035	1.0035	1.0036
23.5	1.0034	1.0035	1.0035	1.0036	1.0036	1.0037
24	1.0035	1.0036	1.0036	1.0037	1.0038	1.0038
24.5	1.0037	1.0037	1.0038	1.0038	1.0039	1.0039
25	1.0038	1.0038	1.0039	1.0039	1.0040	1.0041
25.5	1.0039	1.0040	1.0040	1.0041	1.0041	1.0042
26	1.0040	1.0041	1.0042	1.0042	1.0043	1.0043
26.5	1.0042	1.0042	1.0043	1.0043	1.0044	1.0045
27	1.0043	1.0044	1.0044	1.0045	1.0045	1.0046
27.5	1.0044	1.0045	1.0046	1.0046	1.0047	1.0047
28	1.0046	1.0046	1.0047	1.0048	1.0048	1.0049
28.5	1.0047	1.0048	1.0048	1.0049	1.0050	1.0050
29	1.0049	1.0049	1.0050	1.0050	1.0051	1.0052
29.5	1.0050	1.0051	1.0051	1.0052	1.0052	1.0053
30	1.0052	1.0052	1.0053	1.0053	1.0054	1.0055



10589070




1053950



1131710



1057480 for Grey
1057490 for Orange
1057500 for Yellow
1057510 for Blue
1057520 for Green



2202560 0,5-10 µl
2-20 µl
2202570 5-40 µl
20-200 µl
2202580 40-200 µl
1-5 ml

1-5 ml
200-1000 µl
40-200 µl
5-40 µl
0,5-10 µl



2201860



10589450



1053840



2201520



2202310



1130560



1130510



1131380



1131090



1130550



1053860



1053870



1053860



1054260



1030160



1090610



1053860



1030230



1030020



2x103017



1056350



2201850



1052670



1054030



1056350



2x1030060

2202330



1054730



1053800



1053810



1053810



1053810

SPARE PARTS


Tool
10589070



Cap, red
1058260




2-10 ml
Handle
2203510



Tip Ejector
Part 1
10589600




Piston
Assy
2204890



O-Ring
1033050



Cylinder
1056280



Tip Ejector
Part 2
10589620



Tip Ejector
pusher
1053950



Spring
1131710






Maintenance
Pliers
2900360








Disposable Tip FT-63
Box of 40 pcs, 9402150
Bag of 100 pcs, 9402151




SPARE PARTS 10 ml









8



9



10




11

TIP ORDERING INFORMATION

Code	Finntip	Volume	Qty
9400310	10	0,2 - 10 μ l	1000/bag
9400300	10	0,2 - 10 μ l	10x96/rack
9400130	200 Ext	5 - 200 μ l	10x96/rack
9400260	250 Univ.	0,5 - 250 μ l	10x96/rack
9401250	300	5 - 300 μ l	10x96/rack
9401070	1000	100 - 1000 μ l	200/box
9401110	1000	100-1000 μ l	10x96/rack
9402070	5 ml	1-5 ml	5x54/tray
9402160	10 ml	2-10 ml	5x24/tray



10



200 Ext



250
Universal



300



1000



5 ml



10 ml

TIPS

Product specifications are subject to change without prior notice. Finpipette® and Finntip® are registered trademarks of Thermo Electron Oy.