**Finnpipette® Focus**  
**Single Channel**  
**Variable & Fixed Volume**  
Instructions for Use  
Bedienungsanleitung  
Guide d'utilisation  
Instruccions de uso  
取扱説明書

Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION

Product specifications are subject to change without prior notice. Finnpiquette® and Finntip® are registered trademarks of Thermo Electron Oy.

See the latest version at **[www.thermo.com/finnpiquette](http://www.thermo.com/finnpiquette)**

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**CONTENTS**

PRODUCT DESCRIPTION	4
PIPETTE OPERATION	5
PIPETTING TECHNIQUES	5
CALIBRATION	7
MAINTENANCE	9
TROUBLE SHOOTING	11
PACKAGE	11
SPARE PARTS	48-54

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## Product description

The Finnpiquette FOCUS is an autoclavable pipette. It operates on the air displacement principle and uses detachable, disposable tips.

The adjusted delivery volume is displayed on a numeric readout window on the handle.

The thirteen different models of Finnpiquette FOCUS pipettes with short tip cones cover a volume range from 0.3 µl to 10 ml.

The eleven different models of Finnpiquette FOCUS pipettes with long tip cones cover a volume range from 0.3 µl to 1000 µl.

<b>Order No.</b>	<b>Short</b>	<b>Long</b>	<b>Volume Range</b>	<b>Finntip</b>
4600100	4600150	0,3 µl	to 3 µl	10
4600110	4600160	0,5 µl	to 5 µl	10
4600080	4600170	1 µl	to 10 µl	10
4600010	4600180	1 µl	to 10 µl	250 Universal, 200 Ext
4600120	4600190	3 µl	to 30 µl	50
4600130	4600200	3 µl	to 30 µl	250 Universal, 200 Ext
4600090	4600210	5 µl	to 50 µl	50
4600020	4600220	5 µl	to 50 µl	250 Universal, 200 Ext
4600030	4600230	10 µl	to 100 µl	250 Universal, 200 Ext
4600040	4600240	30 µl	to 300 µl	300
4600050	4600250	100 µl	to 1000 µl	1000
4600060	-	0,5 ml	to 5 ml	5 ml
4600070	-	1 ml	to 10 ml	10 ml

The sixteen different models of Finnpiquette FOCUS Fixed Volume pipettes with short tip cones cover a volume range from 1 µl to 10 ml.

The twelve different models of Finnpiquette FOCUS Fixed Volume pipettes with long tip cones cover a volume range from 1 µl to 1000 µl.

<b>Order No.</b>	<b>Short</b>	<b>Long</b>	<b>Volume Range</b>	<b>Finntip</b>
4601000	4601160	1 µl		10
4601010	4601170	2 µl		10
4601020	4601180	5 µl		10
4601030	4601190	10 µl		10
4601040	4601200	20 µl		250 Universal, 300, 200 Ext
4601050	4601210	25 µl		250 Universal, 300, 200 Ext
4601060	4601220	50 µl		250 Universal, 300, 200 Ext
4601070	4601230	100 µl		250 Universal, 300, 200 Ext
4601080	4601240	200 µl		250 Universal, 300, 200 Ext
4601090	4601250	250 µl		250 Universal, 300
4601100	4601260	500 µl		1000
4601110	4601270	1000 µl		1000
4601120	-	2 ml		5 ml
4601130	-	3 ml		5 ml
4601140	-	5 ml		5 ml
4601150	-	10 ml		10 ml

### Numeric display

The adjusted delivery volume is clearly indicated in the numeric display on the handle.

### Raw materials

The Finnpiquette FOCUS is made of mechanically durable and chemically resistant materials which allow repeated autoclaving of the complete pipette at 121°C (252°F).

The user should verify that the pipette and tip is resistant to the liquid which it is intended to be used.

### Description of tips


Finntips are recommended for use with the Finnpiquette FOCUS.

They are made of natural colour polypropylene, generally regarded as the only contamination free material suitable for tips. Finntips are also autoclavable (121°C / 252°F).

## Pipette operation

### Setting the delivery volume

- Rapid adjustment:** Set the delivery volume using the push button on the top of the pipette. To increase the delivery volume, turn the push button counterclockwise. To decrease the delivery volume, turn it clockwise.
- Fine adjustment:** Set the last digit (0-9) of the delivery volume by turning the fine adjustment ring.
- Make sure that the desired delivery volume clicks into place and that the digits are completely visible in the display window.

The adjustable volume is limited to the specific volume range.

Using excessive force to turn the push button outside the range may jam the mechanism and eventually damage the pipette.

### Tip ejection


To help eliminate the risk of contamination, each pipette is fitted with a tip ejector system.

The tip ejector system consists of a soft-touch tip ejector and specially designed gearing mechanism. To release the tip, point the pipette at suitable waste receptacle and press the tip ejector with your thumb.

### Shelf hanger

You can attach the pipette shelf hanger on a counter, pipette stand or anywhere where you want to hang your pipette.

Clean the area where you plan to attach the shelf hanger. Apply two stickers to the underside of the shelf hanger. Press the shelf hanger firmly into place — on a shelf, counter top or pipette stand. To use, hang the grippy finger rest on the shelf hanger.

## Pipetting techniques

### General

Special care shall be taken to avoid liquid being aspirated into the interior working of the pipette.

Push and release the push button slowly at all times, particularly when working with high viscosity liquids. Never allow the push button to snap back.


Make sure that the tip is firmly attached to the tip cone. Check for foreign particles in the tip.

Before you begin your actual pipetting work, fill and empty the tip 3-5 times with the solution that you will be pipetting. Hold the pipette in an upright position while aspirating liquid. The grippy should rest on your index finger. Make sure that the tips, pipette and solution are at the same temperature.

### Forward technique

Fill a clean reagent reservoir with the liquid to be dispensed.

- Depress the push button to the first stop.
- Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm and slowly release the push button. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
- Deliver the liquid by gently depressing the push button to the first stop. After a delay of about one second, continue to depress the push button all the way to the second stop. This action will empty the tip.
- Release the push button to the ready position.


If necessary, change the tip and continue pipetting.

## Reverse technique

The reverse technique is suitable for dispensing liquids that have a high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes. Fill a clean reagent reservoir with the liquid to be dispensed.


1. Depress the push button all the way to the second stop.
2. Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm and slowly release the push button. This action will fill the tip. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the push button to the first stop. Hold the push button at the first stop. Some liquid will remain in the tip and this should not be included in the delivery.
4. The remaining liquid should either be discarded with the tip or pipetted back into the container.

## Repetitive technique

The repetitive technique offers a rapid and simple procedure for repeated delivery of the same volume. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the push button all the way to the second stop.
2. Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm and slowly release the push button. This action will fill the tip. Withdraw the tip from the liquid touching against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the push button to the first stop. Hold the push button at the first stop. Some liquid will remain in the tip and this should not be included in the delivery.
4. Continue pipetting by repeating steps 2 and 3.


## Pipetting whole blood

(determination in blood glucose determination, for example)

Use steps 1 and 2 of the forward technique to fill the tip with blood.

Wipe the tip carefully with a dry clean tissue.

1. Immerse the tip into the reagent and depress the push button to the first stop, making sure the tip is well below the surface.
2. Release the push button slowly to the ready position. This will fill the tip. Keep the tip in the solution.
3. Depress the push button to the first stop and release slowly. Keep repeating this procedure until the interior wall of the tip is clear.
4. Finally, depress the push button all the way to the second stop to completely empty the tip.


## Changing the grip

### Removal of the grip:

1. Raise the lower end of grip (snap fitting) (3)
2. Take both black side covers out of the grip.

### Assembly of the grip:

1. Put the two side covers into place in the grip.
2. Hold the grip and covers together with fingers and slide the upper end of the grip into the handle. (1 and 2)
3. Turn the lower end of the grip into the handle (snap) and note that the side covers are properly in the handle. (3)

## Calibration

All Finnpipettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. The pipettes are constructed to permit re-adjustment for liquids of different temperature and viscosity.

### Device requirements and test conditions

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range	readable graduation
under 10 $\mu\text{l}$	0.00 1 mg
10-100 $\mu\text{l}$	0.01 mg
above 100 $\mu\text{l}$	0.1 mg

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ( $\pm 0.5^\circ\text{C}$ ) temperature of water, pipette and air between  $15^\circ\text{C}$  to  $30^\circ\text{C}$ . The relative humidity must be above 50%. Especially with volumes under 50  $\mu\text{l}$  the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

### Procedure to check calibration

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume. A new tip is first pre-wetted 3-5 times and a series of ten pipettings is done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume.

#### Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the accuracy (A) and precision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the calculated results are within the selected limits, the adjustment of the pipette is correct.

TABLE 1: Maximum permissible errors according ISO8655

Range	Volume $\mu\text{l}$	Accuracy		Precision	
		$\mu\text{l}$	%	s.d. $\mu\text{l}$	cv%
0,3-3 $\mu\text{l}$	3	$\pm 0.125$	$\pm 4.2$	0.075	2.5
	0,3	$\pm 0.125$	$\pm 42$	0.075	25
0,5-5 $\mu\text{l}$	5	$\pm 0.125$	$\pm 2.5$	0.075	1.5
	0,5	$\pm 0.125$	$\pm 25$	0.075	15
1-10 $\mu\text{l}$	10	$\pm 0.120$	$\pm 1.2$	0.080	0.8
	1	$\pm 0.120$	$\pm 12$	0.080	8.0
3-30 $\mu\text{l}$	30	$\pm 0.50$	$\pm 1.7$	0.20	0.7
	3	$\pm 0.50$	$\pm 17$	0.20	6.7
5-50 $\mu\text{l}$	50	$\pm 0.50$	$\pm 1.0$	0.20	0.4
	5	$\pm 0.50$	$\pm 10$	0.20	4.0
10-100 $\mu\text{l}$	100	$\pm 0.80$	$\pm 0.8$	0.30	0.3
	10	$\pm 0.80$	$\pm 8.0$	0.30	3.0
30-300 $\mu\text{l}$	300	$\pm 4.0$	$\pm 1.3$	1.5	0.5
	30	$\pm 4.0$	$\pm 13$	1.5	5.0
100-1000 $\mu\text{l}$	1000	$\pm 8.0$	$\pm 0.8$	3.0	0.3
	100	$\pm 8.0$	$\pm 8.0$	3.0	3.0
0,5-5 ml	5000	$\pm 40.0$	$\pm 0.8$	15.0	0.3
	500	$\pm 40.0$	$\pm 8.0$	15.0	3.0
1-10 ml	10000	$\pm 60.0$	$\pm 0.6$	30.0	0.3
	1000	$\pm 60.0$	$\pm 6.0$	30.0	3.0

Volume $\mu\text{l}$	Accuracy		Precision	
	$\mu\text{l}$	%	s.d. $\mu\text{l}$	cv%
1	$\pm 0.04$	$\pm 4.0$	0.04	4.0
2	$\pm 0.05$	$\pm 2.5$	0.04	2.0
5	$\pm 0.07$	$\pm 1.4$	0.07	1.4
10	$\pm 0.09$	$\pm 0.9$	0.08	0.8
20	$\pm 0.12$	$\pm 0.6$	0.10	0.5
25	$\pm 0.15$	$\pm 0.6$	0.13	0.5
50	$\pm 0.3$	$\pm 0.6$	0.2	0.4
100	$\pm 0.4$	$\pm 0.4$	0.3	0.3
200	$\pm 0.8$	$\pm 0.4$	0.6	0.3
250	$\pm 1.0$	$\pm 0.4$	0.8	0.3
500	$\pm 1.5$	$\pm 0.3$	1.5	0.3
1000	$\pm 3.0$	$\pm 0.3$	0.3	0.3
2 ml	$\pm 6.0$	$\pm 0.3$	4.0	0.2
3 ml	$\pm 9.0$	$\pm 0.3$	6.0	0.2
5 ml	$\pm 15.0$	$\pm 0.3$	10.0	0.2
10 ml	$\pm 30.0$	$\pm 0.3$	20.0	0.2

## Adjustment

Adjustment is done with the service tool.

1. Place the service tool into the openings of the calibration nut at the top of the handle. Keep fine volume adjustment in place.
2. Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.
3. After adjustment check the calibration according to the instructions above.

## Formulas for calculating results

### Conversion of mass to volume

$$V = (w + e) \times Z$$

$V$  = volume ( $\mu\text{l}$ )  
 $w$  = weight (mg)  
 $e$  = evaporation loss (mg)  
 $Z$  = conversion factor for mg/ $\mu\text{l}$  conversion

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6 mg = 0.2 mg/s).

Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2 mg (10 s x 0.2 mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary.

The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032  $\mu\text{l}/\text{mg}$  at 22°C and 95 kPa. See the conversion table on page 47.

### Accuracy (systematic error)

Accuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0$$

$A$  = accuracy  
 $\bar{V}$  = mean volume  
 $V_0$  = nominal volume

Accuracy can be expressed as a relative value:  $A\% = 100\% \times A / V_0$

### Precision (random error)

Precision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$S = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

$s$  = standards deviation  
 $\bar{v}$  = mean volume  
 $n$  = number of measurements

Standard deviation can be expressed as a relative value (CV)  $CV = 100\% \times S / \bar{V}$

## Maintenance

When the Finnpiquette FOCUS is not in use, make sure it is stored in an upright position. We recommend a Finnpiquette stand for this purpose.

The part # refer to exploded views beginning at page 48 and 52 Long model.

### Daily checking

The pipette should be checked at the beginning of each day for dust and dirt on the outside surfaces of the pipette. Particular attention should be paid to the tip cone. No other solvents except 70 % ethanol should be used to clean the pipette.

### Short-term service


If the pipette is used daily it should be checked and lubricated every three months.

#### Pipettes up to 1000 µl

1. Remove the lower tip ejector sleeve 13 by pulling it out (snap fitting).
2. Push down the ejector button and take the upper ejector top 11 between fingers.
3. Release the ejector button to up position, but hold the tip ejector top 11 in down position with the fingers.
4. Remove the tip cone by pulling it out (snap fitting).
5. Pull out the piston 14 and the spring 15 but not spring 16.
6. Clean the piston 14 and the piston spring 15 with a dry napless cloth.
7. Check the tip cone for foreign particles.
8. Grease the cleaned parts with the lubricant that comes with the pipette.
9. Reassemble the pipette components in reverse order.

#### 5 ml and 10 ml pipettes / 2-10 ml pipettes

1. Pull down the tip ejector sleeve.
2. Note that the tip ejector button is in fully up position.
3. Remove the tip ejector by pulling it out.
4. Remove the ejector sleeve by pulling it to two parts (snap fitting).
5. Remove the cylinder 14 by pushing the piston to the down position and push the snap joints to release and remove the cylinder 14.
6. Push the piston forward and clean it with a dry napless cloth.
7. Check the tip cone for foreign particles.
8. Grease the cleaned parts with the lubricant that comes with the pipette.
9. Reassemble the pipette components in reverse order.

### Long-term service

If the pipette is used daily it should be serviced every six months.

#### Pipettes up to 1000µl

1. Remove the lower tip ejector sleeve 13 by pulling it out (snap fitting).
2. Push down the ejector button and take the upper ejector top 11 between fingers.
3. Release the ejector button to up position, but hold the tip ejector top 11 in down position with the fingers.
4. Remove the tip cone by pulling it out (snap fitting).
5. Pull out the piston 14 and the spring 15.
6. Remove spring hatch 25 by turning it. Take out the spring 16.
7. Remove the rest of the parts from tip cone with the piston.
8. Clean the parts with a dry napless cloth.
9. Check the tip cone for foreign particles.
10. Grease the cleaned parts with the lubricant that comes with the pipette.

#### 5 ml and 10 ml pipettes / 2-10 ml pipettes

The long-term service for 5ml and 10ml pipettes is the same as short-term service.

## Reassembly:

### 1-10 $\mu\text{l}$ , 0.5-5 $\mu\text{l}$ and 0.3-3 $\mu\text{l}$ pipettes:

1. Put the spring hatch 25, the spring 16, spring support 17, (part 26 only long model) and the tube 18 back on the piston.
2. Slide the bigger O-ring 19, smaller O-ring 20 and tube 21 on the piston.
3. Slide the small spring 22, spring support 23 and O-ring 24 on the tube 21.
4. Carefully slide the entire assembly into the tip, then push the spring 16 until it snaps.
5. Pull out the piston, insert the spring 15 on piston, and put the piston back into the tip cone.
6. Take the upper ejector top 11 between fingers and pull it down.
7. Push the tip cone into the handle while holding the ejector top 11 down.
8. Reassemble the lower tip ejector sleeve 13 (snap fitting).

### 5-50 $\mu\text{l}$ and 3-30 $\mu\text{l}$ pipettes / 20-50 $\mu\text{l}$ pipettes:

1. Put the spring hatch 25, the spring 16, spring support 17 and the tube 18 back on the piston.
2. Slide the bigger O-ring 19, smaller O-ring 20 and O-ring support 21 on the piston.
3. Slide the small spring 22 on the piston.
4. Carefully slide the entire assembly into the tip, then push the spring 16 until it snaps.
5. Pull out the piston, insert the spring 15 on piston, and put the piston back into the tip cone.
6. Take the upper ejector top 11 between fingers and pull it down.
7. Push the tip cone into the handle while holding the ejector top 11 down.
8. Reassemble the lower tip ejector sleeve 13 (snap fitting).

### 10-100 $\mu\text{l}$ pipette / 100 $\mu\text{l}$ pipettes:

1. Put the spring hatch 25, the spring 16, spring support 17 and the O-ring 20 back on the piston.
2. Carefully slide the entire assembly into the tip, then push the spring 16 until it snaps.
3. Pull out the piston, insert the spring 15 on piston, and put the piston back into the tip cone.
4. Take the upper ejector top 11 between fingers and pull it down.
5. Push the tip cone into the handle while holding the ejector top 11 down.
6. Reassemble the lower tip ejector sleeve 13 (snap fitting).

### 30-300 $\mu\text{l}$ pipette / 200-300 $\mu\text{l}$ pipettes:

1. Put the spring 16, spring support 17 and the O-ring 20 back on the piston.
2. Carefully slide the entire assembly into the tip, then push the spring 16 until it snaps.
3. Pull out the piston, insert the spring 15 on piston, and put the piston back into the tip cone.
4. Take the upper ejector top 11 between fingers and pull it down.
5. Push the tip cone into the handle while holding the ejector top 11 down.
6. Reassemble the lower tip ejector sleeve 13 (snap fitting).

### 100-1000 $\mu\text{l}$ pipette / 500-1000 $\mu\text{l}$ pipettes:

1. Put the spring 16, spring support 17 and the O-ring 20 back on the piston.
2. Put the spring 15 on the piston.
3. Carefully slide the entire assembly into the tip cone.
4. Push the spring 16 until it snaps.
5. Take the upper ejector top 11 between fingers and pull it down.
6. Push the tip cone into the handle while holding the ejector top 11 down.
7. Reassemble the lower tip ejector sleeve 13 (snap fitting).

## Sterilization

The entire pipette can be sterilized by autoclaving it at 121°C (252°F) (minimum 20 minutes). No special preparations are needed for autoclaving. You can use steam sterilization bags if needed.

After autoclaving the pipette must be cooled to room temperature for at least two hours. Before pipetting, make sure that the pipette is dry. We recommend that you check the calibration after every sterilization cycle to achieve the best possible accuracy.

## Trouble shooting

The table below lists possible problems and their solutions.

<b>Defect</b>	<b>Possible reason</b>	<b>Solution</b>
<i>Leakage</i>	<i>Tip incorrectly attached</i>	<i>Attach firmly</i>
	<i>Foreign particles between tip and tip cone</i>	<i>Clean tip cones attach new tips</i>
	<i>Foreign particles between the piston, the O-ring and the cylinder</i>	<i>Clean and grease O-ring and cylinder.</i>
	<i>Insufficient amount of grease on cylinder and O-ring</i>	<i>Grease accordingly</i>
	<i>O-ring damaged</i>	<i>Change the O-ring</i>
<i>Inaccurate dispensing</i>	<i>Incorrect operation</i>	<i>Follow instructions carefully</i>
	<i>Tip incorrectly attached</i>	<i>Attach firmly</i>
	<i>Calibration altered: caused by misuse, for example</i>	<i>Recalibrate according to instructions</i>
<i>Inaccurate dispensing with certain liquids</i>	<i>Unsuitable calibration</i>	<i>Recalibrate with the liquids in question</i>
	<i>High viscosity liquids may require recalibration</i>	

## Package

The Finnpiquette FOCUS is shipped in a specially designed package containing the following items:

1. The Finnpiquette
2. Service tool
3. Finntip sample
4. Tube of grease
5. Instruction manual
6. Calibration certificate
7. Shelf hanger

### **CAUTION!**

***The Finnpiquette is designed to allow easy in-lab service. If you would prefer to have us or your local representative service your pipette, please make sure that the pipette has been decontaminated before you send it to us.***

***Please note that the postal authorities in your country may prohibit or restrict the shipment of contaminated material by mail.***

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Thermo Electron Corporation