**Fisherbrand**<sup>®</sup> *Elite*

## Single Channel & Multichannel pipette

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## English

### CONTENTS

PRODUCT DESCRIPTION	3
PIPETTE OPERATION	3
PIPETTING TECHNIQUES	4
CALIBRATION AND ADJUSTMENT	4
MAINTENANCE	7
TROUBLE SHOOTING	9
PACKAGE	9
SPARE PARTS	27-34

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## Product description

The **Fisherbrand**<sup>®</sup> Elite is a continuously adjustable, general purpose micropipette for sampling and dispensing accurate volumes of liquid.

It operates on an air displacement principle (i.e. an air interface) and uses detachable, disposable tips.

The adjusted delivery volume is displayed digitally on a readout window in the handle.

The ten different models of **Fisherbrand**<sup>®</sup> Elite pipettes cover a volume range from 0,2 µl to 10 ml.

Order No.	Volume Range	SureOne
FB77100	0,2 µl to 2 µl	<b>10</b>
FB77102	0,5 µl to 5 µl	<b>10</b>
FB77104	1 µl to 10 µl	<b>10</b>
FB77106	2 µl to 20 µl	<b>200</b>
FB77108	5 µl to 50 µl	<b>200</b>
FB77110	10 µl to 100 µl	<b>200</b>
FB77112	20 µl to 200 µl	<b>200</b>
FB77114	100 µl to 1000 µl	<b>1250</b>
FB77116	0,5 ml to 5 ml	<b>5 ml</b>
FB77118	1 ml to 10 ml	<b>10 ml</b>

The eight different models of **Fisherbrand**<sup>®</sup> Elite Multichannel pipettes cover a volume range from 1 µl to 300 µl.

Order No.	Channel	Volume Range	SureOne
FB77120	8	1 µl to 10 µl	<b>10</b>
FB77124	8	5 µl to 50 µl	<b>200</b>
FB77122	8	10 µl to 100 µl	<b>200</b>
FB77126	8	30 µl to 300 µl	<b>300</b>
FB77128	12	1 µl to 10 µl	<b>10</b>
FB77132	12	5 µl to 50 µl	<b>200</b>
FB77130	12	10 µl to 100 µl	<b>200</b>
FB77134	12	30 µl to 300 µl	<b>300</b>

### Digital display

The adjusted delivery volume is clearly indicated in the large digital display on the handle.

### Raw materials

The **Fisherbrand**<sup>®</sup> Elite is made of mechanically durable and chemically resistant materials.

### Description of tips

SureOne pipette tips are recommended for use with the **Fisherbrand**<sup>®</sup> Elite.

They are made of virgin natural colour polypropylene, generally regarded as the only contamination free material suitable for tips. SureOne standard tips are autoclavable (121°C), the filter tips are not autoclavable.

## Pipette operation

### Setting the delivery volume

1. Set the delivery volume using the push button on the top of the pipette. To increase the delivery volume, turn the push button counterclockwise. To decrease the delivery volume, turn it clockwise.
2. Make sure that the desired delivery volume clicks into place.
3. Do not set volumes outside the pipette's specified volume range. Using excessive force to turn the push button outside the range may jam the mechanism and eventually damage the pipette.

## Tip ejection

To help eliminate the risk of contamination, each pipette is fitted with a tip ejector system. To release the tip, point the pipette at suitable waste receptacle and press the tip ejector with your thumb.

## Pipetting techniques

Push and release the push button slowly at all times particularly when working with high viscosity liquids. Never allow the push button to snap back.

Make sure that the tip is firmly attached to the tip cone. Check for foreign particles in the tip.

Before you begin your actual pipetting work, fill and empty the tip 2-3 times with the solution that you will be pipetting. Hold the pipette in an upright position while aspirating liquid. The pipette should be positioned in your hand with the finger rest against your index finger. Make sure that the tips, pipette and solution are at the same temperature.

### Forward technique

Fill a clean reagent reservoir with the liquid to be dispensed.

- 1 Depress the push button to the first stop.
- 2 Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm and slowly release the push button. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
- 3 Deliver the liquid by gently depressing the push button to the first stop. After a delay of about one second, continue to depress the push button all the way to the second stop. This action will empty the tip.
- 4 Release the push button to the ready position. If necessary, change the tip and continue pipetting.

### Reverse technique

The reverse technique is suitable for dispensing liquids that have a high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the push button all the way to the second stop.
2. Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm, and slowly release the push button. This action will fill the tip. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the push button to the first stop. Hold the push button at the first stop. Some liquid will remain in the tip and this should not be included in the delivery.
4. The remaining liquid should either be discarded with the tip or dispensed back into the reservoir.

## Calibration and adjustment

All **Fisherbrand**® Elite pipettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water using the forward pipetting technique. It should be noted that the use of other pipetting techniques may affect the calibration results. The pipettes are constructed to permit re-adjustment for other pipetting techniques or liquids of different temperature and viscosity.

## Device requirements and test conditions

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range	readable graduation
under 10 $\mu\text{l}$	0.00 1 mg
10-100 $\mu\text{l}$	0.01 mg
above 100 $\mu\text{l}$	0.1 mg

Test liquid: Water, distilled or deionized, "grade 3" water conforming to ISO 3696. Tests are done in a draft-free room at a constant ( $\pm 0.5^\circ\text{C}$ ) temperature of water, pipette and air between  $15^\circ\text{C}$  to  $30^\circ\text{C}$ . The relative humidity must be above 50%. Especially with volumes under 50  $\mu\text{l}$  the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

## Procedure to check calibration

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume. A new tip is first pre-wet 3-5 times and a series of ten pipettings is done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume. Use of forward pipetting technique is recommended. The maximum permissible errors are designed for forward method.

### Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the inaccuracy (A) and imprecision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the calculated results are within the selected limits, the adjustment of the pipette is correct.

TABLE 1: Maximum permissible errors according ISO8655

Range	Volume $\mu\text{l}$	Inaccuracy		Imprecision	
		$\mu\text{l}$	%	s.d. $\mu\text{l}$	cv%
0,2-2 $\mu\text{l}$	2	$\pm 0.080$	$\pm 4$	0.040	2.0
	0.2	$\pm 0.080$	$\pm 40$	0.040	20.0
0,5-5 $\mu\text{l}$	5	$\pm 0.125$	$\pm 2.5$	0.075	1.5
	0.5	$\pm 0.125$	$\pm 25$	0.075	15
1-10 $\mu\text{l}$	10	$\pm 0.120$	$\pm 1.2$	0.080	0.8
	1	$\pm 0.120$	$\pm 12$	0.080	8.0
2-20 $\mu\text{l}$	20	$\pm 0.20$	$\pm 1.0$	0.10	0.5
	2	$\pm 0.20$	$\pm 10.0$	0.10	5.0
5-50 $\mu\text{l}$	50	$\pm 0.50$	$\pm 1.0$	0.20	0.4
	5	$\pm 0.50$	$\pm 10$	0.20	4.0
10-100 $\mu\text{l}$	100	$\pm 0.80$	$\pm 0.8$	0.30	0.3
	10	$\pm 0.80$	$\pm 8.0$	0.30	3.0
20-200 $\mu\text{l}$	200	$\pm 1.60$	$\pm 0.8$	0.60	0.3
	20	$\pm 1.60$	$\pm 8.0$	0.60	3.0
100-1000 $\mu\text{l}$	1000	$\pm 8.0$	$\pm 0.8$	3.0	0.3
	100	$\pm 8.0$	$\pm 8.0$	3.0	3.0
0,5-5 ml	5000	$\pm 40.0$	$\pm 0.8$	15.0	0.3
	500	$\pm 40.0$	$\pm 8.0$	15.0	3.0
1-10 ml	10000	$\pm 60.0$	$\pm 0.6$	30.0	0.3
	1000	$\pm 60.0$	$\pm 6.0$	30.0	3.0

Range	Channel	Volume $\mu\text{l}$	Inaccuracy		Imprecision	
			$\mu\text{l}$	%	s.d. $\mu\text{l}$	cv%
1-10 $\mu\text{l}$	8, 12	10	$\pm 0.24$	$\pm 2.4$	0.16	1.6
		1	$\pm 0.24$	$\pm 24$	0.16	16
5-50 $\mu\text{l}$	8, 12	50	$\pm 1.0$	$\pm 2.0$	0.4	0.8
		5	$\pm 1.0$	$\pm 20$	0.4	8.0
10-100 $\mu\text{l}$	8, 12	100	$\pm 0.80$	$\pm 0.8$	0.30	0.3
		10	$\pm 0.80$	$\pm 8.0$	0.30	3.0
30-300 $\mu\text{l}$	8, 12	300	$\pm 8.0$	$\pm 2.7$	3.0	1.0
		30	$\pm 8.0$	$\pm 26.7$	3.0	10.0

## Adjustment

Adjustment is done with the service tool.

1. Place the service tool into the openings of the calibration nut at the top of the handle.
2. Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.
3. After adjustment check the calibration according to the instructions above.

## Formulas for calculating results

### Conversion of mass to volume

$$V = (w + e) \times Z$$

V = volume ( $\mu\text{l}$ )

w = weight (mg)

e = evaporation loss (mg)

Z = conversion factor for  $\mu\text{l}/\text{mg}$  conversion

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6 mg = 0.2 mg/s).

Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2 mg (10 s x 0.2 mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary.

The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032  $\mu\text{l}/\text{mg}$  at 22°C and 95 kPa. See the conversion table on page 26.

### Inaccuracy (systematic error)

Inaccuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0 \quad \begin{array}{l} \frac{A}{\bar{V}} = \text{inaccuracy} \\ \bar{V} = \text{mean volume} \\ V_0 = \text{nominal volume} \end{array}$$

Inaccuracy can be expressed as a relative value:  $A\% = 100\% \times A / V_0$

### Imprecision (random error)

Imprecision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$S = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}} \quad \begin{array}{l} s = \text{standards deviation} \\ \bar{V} = \text{mean volume} \\ n = \text{number of measurements} \end{array}$$

Standard deviation can be expressed as a relative value (CV)

$$CV = 100\% \times S / \bar{V}$$

# Maintenance

When the **Fisherbrand**® Elite is not in use, make sure it is stored in an upright position. We recommend a pipette stand for this purpose.

For a complete part number listing refer to the pipette exploded views beginning on page 27.

## Short-term service

The pipette should be checked at the beginning of each day for dust and dirt on the outside surfaces of the pipette.

Particular attention should be paid to the tip cone. No other solvents except 70 % ethanol should be used to clean the pipette.

## Long-term service, single channel pipettes

If the pipette is used daily it should be checked every three months.

The servicing procedure starts with the disassembly of the pipette.

### 1-1000 µl pipettes

1. Press the tip ejector.
2. Rotate the tip ejector 11 counterclockwise and pull it out.
3. Turn the tip cone counterclockwise with the service tool to remove.
4. Pull out the piston and other parts. Push out with piston the rest of the piston assembly. Then turn the tip cone upside down and tap all parts from tip cone. Remember keep all parts in order on table for reassembly.
5. Clean the piston, the piston spring and the O-rings with a dry napless cloth.
6. Check the tip cone for foreign particles.
7. Grease the cleaned parts with the lubricant that comes with the pipette.
8. Reassemble the pipette components.


**0.2-2µl, 0.5-5µl & 1-10 µl:** First, slide spring 22, o-ring support 23 and o-ring 24 on the tube 21. With the 0.2-2µl model insert the tube 27 into the tube 21. Then slide the spring 13, spring support 16 and tubes 17 and 18, bigger o-ring 19 and smaller o-ring 20 back on the piston. Compress the spring with fingers by pressing piston and spring support 16 against each other and slide the tube 21 with rest of the parts on the piston. Hold the spring compressed and carefully slide the entire assembly into the tip cone and release the spring.

**2-20 µl & 5-50 µl:** Slide the spring 13, spring support 16 and tubes 17 and 18, bigger o-ring 19 and smaller o-ring 20 back on the piston.

Compress the spring with fingers by pressing piston and spring support 16 against each other and slide the bigger o-ring 19, smaller o-ring 20, spring support 21 and the spring 22 (smaller diameter against spring support 21) on the piston. Hold the spring compressed and carefully slide the entire assembly into the tip cone and release the spring.

**10-100 µl & 20-200 µl:** Slide the spring 13, spring support 16 and o-ring 17 back on the piston. Slide the entire assembly into the tip cone.


**100-1000µl:** Put the o-ring 17 and support ring 16 to the tip cone. Slide the spring 13 on the piston and slide the entire assembly into the tip cone.

- All:** Put the spring 15 and support 14 on top of the tip cone and carefully insert the tip cone assembly to the handle and turn it tight by hand.
- Reassemble the tip ejector.

### 0.5-5ml & 1-10 ml pipettes

- Press the tip ejector.
- Rotate the tip ejector 10 counterclockwise to open it.
- Disassemble the lower part of the tip ejector 14 (snap fitting).
- Turn the cylinder 13 counterclockwise and pull out the tip cone assembly.
- Remove the cylinder 13 by pressing the snaps fittings of the cylinder.
- Clean and regrease the O-ring 12 and cylinder 13.
- Assemble the parts in the opposite order of disassembly.

### Long-term service, multichannel pipettes

If the pipette is used daily it should be checked and lubricated every three months.

- Place the service tool head#1 between ring 15 and ejector 23. Push the tool until the parts snap from each other.
- Check that the ejector lever is in the up position and pull down the tip ejector part of the module.  
Place the service tool head#2 in the hole of adapter tube 46.
- Open the upper end of the tip ejector slightly and remove the tip ejector.
- Unscrew the module from the handle.
- Pull out the module spring 19 and clip 22.
- Press the spring 13 and remove the locking pieces 12 from the groove. Remove the spring 13.
- Take off the locking claws 44 and 45 and pull out the adapter tube and tube 43.
- Use a screwdriver to remove the four screws in the module cover and lift off the cover.
- Remove the piston bar and clean the pistons and tip cones with a dry napless cloth.
- If needed, service the tip cones:

**30-300  $\mu$ l, 10-100 $\mu$ l & 5-50  $\mu$ l:** Open the tip cone by carefully releasing the cover ring from its snap joint with the screwdriver. Remove all the parts from the tip cone. Clean all the parts. If needed, replace the o-rings. Take one piston. Slide the spring 33, cover ring 32 (larger hole), spring 34, support ring 35, (o-ring 37 bigger 5-50 $\mu$ l) and o-ring 36 (smaller) onto the piston. Grease the o-ring with the lubricant provided in the pipette package. Slide all the parts into the tip cone and close the snap joint of the cover ring.

**1-10  $\mu$ l:** Open the tip cone by carefully releasing the cover ring from its snap joint with the screwdriver. Remove all the parts from the tip cone. Clean all the parts. If needed, replace the o-rings. Take one piston. Slide spring 33, cover ring 32 (larger hole), support 35, o-ring 36 (bigger), o-ring 37 (smaller) and o-ring support 38 onto the piston. Then slide spring 39, spring support 40 (sharp edges first) and o-ring 41 onto the o-ring support 38. Grease the o-rings with the lubricant provided in the pipette package. Slide all the parts into the tip cone and close the snap joint of the cover ring.

- Install the piston bar with pistons and tip cones in the cover and close the cover with the four screws. Insert the clip 22.
- Place the adapter tube and tube 43 on the neck of the module and insert the locking claws 44 and 45. Insert the module spring 19.
- Insert spring 13 and locking pieces 12 to the piston rod 16.
- Place the tip ejector on the module. Push the spring 19 inside tip ejector parts and close the upper end of the ejector and keep closed with fingers.
- Screw the module into the handle and tighten with service tool head#2.
- Push the tip ejector lever down, until you hear a "click".

## Service Instructions for Multichannel Pipette Tip Cones

To ensure performance between all channels in a multichannel pipette you must replace all channels at the same time. Do not replace just the channels that are in need of replacement. Don't mix tip cones of different packages, because one bag contains a matched set of tip cones.

### Sterilization

The entire pipette can be sterilized by autoclaving it at 121°C (252°F) (2 ata) (minimum 20 minutes). No special preparations are needed for autoclaving. You can use steam sterilization bags if needed.

After autoclaving the pipette must be cooled to room temperature for at least two hours. Before pipetting, make sure that the pipette is dry. We recommend that you check the calibration after every sterilization cycle to achieve the best possible accuracy.

## Trouble shooting

The table below lists possible problems and their solutions.

<b>Defect</b>	<b>Possible reason</b>	<b>Solution</b>
<i>Leakage</i>	<i>Tip incorrectly attached</i>	<i>Attach firmly</i>
	<i>Foreign particles between tip and tip cone</i>	<i>Clean tip cones attach new tips</i>
	<i>Foreign particles between the piston, the O-ring and the cylinder</i>	<i>Clean and grease O-ring and cylinder.</i>
	<i>Insufficient amount of grease on cylinder and O-ring</i>	<i>Grease accordingly</i>
<i>Inaccurate dispensing</i>	<i>O-ring damaged</i>	<i>Change the O-ring</i>
	<i>Incorrect operation</i>	<i>Follow instructions carefully</i>
	<i>Tip incorrectly attached</i>	<i>Attach firmly</i>
<i>Inaccurate dispensing with certain liquids</i>	<i>Calibration altered: caused by misuse, for example</i>	<i>Recalibrate according to instructions</i>
	<i>Unsuitable calibration</i>	<i>Recalibrate with the liquids in question</i>
	<i>High viscosity liquids may require recalibration</i>	

## Package

The **Fisherbrand®** Elite is shipped in a specially designed package containing the following items:

1. The **Fisherbrand®** Elite pipette
2. Service tool
3. Multichannel service tool
4. Tube of grease (**Order No. 2203130**)
5. Instruction manual
6. Calibration certificate

### CAUTION!

**The Fisherbrand® Elite is designed to allow easy in-lab service. If you would prefer to have us or your local representative service your pipette, please make sure that the pipette has been decontaminated before you send it to us.**

**Please note that the postal authorities in your country may prohibit or restrict the shipment of contaminated material by mail.**

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Fisher Scientific GTF AB
Fisher Scientific UK Ltd.